



Runs on easy replaceable rechargeable / non rechargeable AA batteries (Batteries not supplied with unit)

Padded headband provides comfortable head grip for large wearing

Small, medium & large Diopter

on/off Switch

Diopter adjustment knob for Small, Medium & Large Pupil

Illumination control knob

illumination Warm white LED

Head adjustable Knob

Sealed optics

Knob To Hold The Teaching Mirror

Mirror tilt control

Four integrated filters
 Green, Blue
 Orange, White

IPD adjustment
 49 to 75 mm

Accessories



20D Lens



Scleral Depressor (Set of 3)



Teaching Mirror



Carrying Case

